



Rat beta 2-microglobulin reference serum (DAGA-686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	β 2 microglobulin also known as B2M is a component of MHC class I molecules, which are present on all nucleated cells (excludes red blood cells). In humans, the β 2 microglobulin protein is encoded by the B2M gene.
Specificity	Rat
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Rat;Beta 2-Microglobulin;Serum
----------	--------------------------------