



Rat prealbumin reference serum (DAGA-700)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Prealbumin (PAB) is a maker of nutritional status and inflammation of potential prognostic value in acute heart failure (HF).PAB can be used as an indicator of liver damage caused by asphyxia in neonates, and can be used to assess the degree of asphyxia
Specificity	Rat
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords	Rat;Prealbumin;Serum
-----------------	----------------------