



Rat vitamin D binding protein reference serum (DAGA-705)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Vitamin D-binding protein, also/originally known as gc-globulin (group-specific component), is a protein that in humans is encoded by the GC gene. Vitamin D-binding protein belongs to the albumin gene family, together with human serum albumin and alpha-fetoprotein. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues.
Specificity	Rat
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 ml
Preservative	0.1% Sodium Azide
Storage	Frozen -20°C

BACKGROUND

Keywords Rat; Vitamin D Binding Protein; Serum