



CA 125 antigen (Control grade,>98%) (DAGA-751)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CA-125 has found application as a tumor marker or biomarker that may be elevated in the blood of some patients with specific types of cancers, or other benign conditions. MUC16 (CA125) is a serum marker that is used routinely in gynecologic practice to monitor patients with ovarian cancer.
Purity	>98% by SDS-PAGE
Conjugate	Unconjugated
Applications	Control
Format	liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS pH 7.2
Preservative	None
Storage	Store at 4°C for a maximum of 12 days. Store at -20°C for long term storage. Avoid freeze/thaw cycles

BACKGROUND

Keywords	MUC16;mucin 16;cell surface associated;CA125;mucin-16;CA125 ovarian cancer antigen;ovarian carcinoma antigen CA125;ovarian cancer-related tumor marker CA125
-----------------	--