



CA 72-4 antigen (DAGA-770)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CA 72-4 is a tumor marker primarily associated with gastric and ovarian cancer, it is also associated with esophageal cancer, pancreatic carcinoma and other tumors.
Purity	>50–80 % by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Concentration	50–500 kU/mL
Size	100kU
Buffer	PBS, pH7.2–7.6
Preservative	0.05% Sodium Azide
Storage	Store at 2–8°C short term and at -20°C long term. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	CA 72-4;CA72-4;Cancer Antigen 72-4 protein;CA72-4 protein
----------	---