



Recombinant PGI antigen (>90%) (DAGA-784)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pepsinogen I is the pro-form of human pepsin A-5. It is produced in the stomach and secreted into the gastric lumen, but a small amount also into the blood. Alterations in the serum pepsinogen concentration may indicate gastric abnormalities
Purity	>90% by SDS-PAGE
Conjugate	TBD
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Buffer	Saline, pH 7.4, containing 3 % sucrose
Preservative	None
Storage	Store at 2-8°C

BACKGROUND

Keywords	PGA3;pepsinogen 3;group I (pepsinogen A);pepsin A-3;pepsin A;pepsinogen-3;pepsinogen A3;Pepsinogen I;pg1
-----------------	--