



Recombinant PGI antigen (>90%) (DAGA-785)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Pepsinogen II is the pro-form of human gastricsin. It is produced in the stomach and secreted into the gastric lumen, but a small amount also into the blood. Alterations in the serum pepsinogen concentration may indicate gastric abnormalities.
Purity	>90 % by SDS-PAGE
Conjugate	TBD
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Buffer	Saline, pH 7.4, containing 3 % sucrose
Preservative	None
Storage	Store at 2-8°C

BACKGROUND

Keywords	PGC;progastricsin (pepsinogen C);PEPC;PGII;gastricsin;pepsin C;pepsinogen C;preprogastricsin;pepsinogen group II;Pg2
-----------------	--