



## PSA protein ( $\geq 96\%$ ) (DAGA-790)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Prostate-specific antigen (PSA) is a protein produced by the cells of the prostate gland. PSA is present in small quantities in the serum of men with healthy prostates, but is often elevated in the presence of prostate cancer and in other prostate disorders.
<b>Purity</b>	$\geq 96\%$ by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA; WB
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	2-8°C

### BACKGROUND

<b>Keywords</b>	Alternative namesAPS;Gamma seminoprotein;Gamma-seminoprotein;hK3;Kallikrein 3;Kallikrein related peptidase 3;Kallikrein-3;KLK 3;KLK2A1;KLK3;KLK3_HUMAN;P 30 antigen;P-30 antigen;P30 antigen;Prostate Specific Antigen;Prostate-specific antigen
-----------------	--