



## CA 242 from Metastatic liver (Immunogen grade) (DAGA-800)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	CA242 antigen is defined by a monoclonal antibody, which was obtained by immunizing mice with a human colorectal carcinoma cell line. A high positive correlation between CA242 values and serum bilirubin or alkaline phosphatase levels has also been earlier reported.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Immunogen, Controls
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	10 kU
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.2 – 7.4 with 3% Sucrose
<b>Preservative</b>	0.05% Sodium Azide
<b>Storage</b>	≤-20°C. Avoid repeated freezing and thawing.

### BACKGROUND

**Keywords** Cancer Antigen 242;CA242;Liver;Tumor Markers