



Recombinant Human HER2 [His] (DAGA-813)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	70.1 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM Citrate, 8% Sucrose, 0.05% Tween 80, pH 4.5
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	HER2; Receptor tyrosine-protein kinase erbB-2; ERBB2
----------	--

GENE INFORMATION

UniProt ID	P04626
------------	--------