



Human Thyroglobulin Biotinylated (DAGA-826B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 80% as determined by SDS-PAGE
Conjugate	Biotin
Applications	N/A
Format	Liquid
Size	100 µg, 1 mg
Buffer	20mM HEPES buffer pH7.6, 150mM NaCl and 40% Sucrose (w/v)
Preservative	None
Storage	Store at 2°C to 8°C short term or frozen below -20°C

BACKGROUND

Keywords	Thyroglobulin; Tg; AITD3
----------	--------------------------