



Myeloperoxidase antigen isoform C (>96%) (DAGA-876)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Myeloperoxidase(MPO) is a peroxidase enzyme (EC 1.11.1.7) most abundantly present in neutrophil granulocytes (a subtype of white blood cells). It is a lysosomal protein stored in azurophilic granules of the neutrophil. MPO has a heme pigment, which causes it
Purity	>96% by SDS-PAGE
Conjugate	Unconjugated
Applications	Controls, Calibrators & Standards Immunoassays Antibody Production Clinical Chemistry Testing/Assay Validation Life Science ELISA, Blotting, Etc.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg, 1 mg
Buffer	50 mM Sodium Acetate, pH 6.0, 100 mM Sodium Chloride
Preservative	None
Storage	2-8°C

BACKGROUND

Keywords	MPO;AFBP;IBP1;PP12;IGF-BP25;hIGFBP-1;IGFBP1;insulin-like growth factor-binding protein 1;IBP-1;IGFBP-1;OTTHUMP00000159561;OTTHUMP00000207990;OTTHUMP00000214348;binding protein-25;binding protein-26;binding protein-28;placental protein 12
-----------------	---