



Human cardiac troponin I subunit antigen (>95%) (DAGA-884)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Human troponin I is presented in cardiac muscle tissue by a single isoform with molecular weight 23876 Da and it consists of 209 amino acid residues. cTnI molecule contains two serines in the 22 and 23 positions. Both amino acid residues can be phosphorylated in vivo by protein kinase A, so four forms of protein-one dephospho, two monophospho and one bisphospho-can coexist in the cell.
Purity	>95% by SDS-PAGE
Conjugate	Unconjugated
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	75 mM Tris, pH 8.0, 15 mM mercaptoethanol, 10mM EDTA, 6 M urea
Preservative	None
Storage	< -20°C or colder long term

BACKGROUND

Keywords Troponin I;Cardiac Troponin 1;cTnI