



# Human apolipoprotein AI/B free serum (DAGA-953)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Apolipoprotein B is the primary apolipoprotein of chylomicrons, VLDL, IDL, and LDL particles (LDL - known commonly by the misnomer "bad cholesterol" when in reference to both heart disease and vascular disease in general), which is responsible for carrying fat molecules (lipids), including cholesterol, around the body (within the water outside cells) to all cells within all tissues.
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Specific methodologies have not been tested using this product.
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Short term store at 2–8°C. Long term store at -20°C. Avoid multiple freeze/thaw cycles.

## BACKGROUND

<b>Keywords</b>	Apolipoprotein B;Apolipoprotein antigen;Apolipoprotein protein
-----------------	--