



Human apolipoprotein AII (>95%) (DAGA-954)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Apolipoprotein B is the primary apolipoprotein of chylomicrons, VLDL, IDL, and LDL particles (LDL - known commonly by the misnomer "bad cholesterol" when in reference to both heart disease and vascular disease in general), which is responsible for carrying fat molecules (lipids), including cholesterol, around the body (within the water outside cells) to all cells within all tissues.
Purity	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	Specific methodologies have not been tested using this product.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	10 mM Ammonium Bicarbonate, pH 7.4
Preservative	None
Storage	Store at < -70°C. Keep concentration > 1mg/mL as precipitation may occur below 1 mg/mL. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Apolipoprotein B;Apolipoprotein antigen;Apolipoprotein protein