



## Human apolipoprotein C1 (>95%) (DAGA-957)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Antigen Description</b>	Apolipoprotein A1 is a protein that in humans is encoded by the APOA1 gene. It has a specific role in lipid metabolism. The text in a recent report suggested that APOA3 mRNA is regulated by endogenously expressed antisense RNA.
<b>Purity</b>	>95% by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Specific methodologies have not been tested using this product
<b>Format</b>	Lyophilized
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Buffer</b>	10 mM NH <sub>4</sub> HCO <sub>3</sub> , pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C. Avoid multiple freeze/thaw cycles.

### BACKGROUND

**Keywords** Apolipoprotein A3;Apolipoprotein antigen;Apolipoprotein protein