



Human apolipoprotein CII (>95%) (DAGA-958)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Apolipoprotein A1 is a protein that in humans is encoded by the APOA1 gene. It has a specific role in lipid metabolism. The text in a recent report suggested that APOA4 mRNA is regulated by endogenously expressed antisense RNA.
Purity	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	Specific methodologies have not been tested using this product
Format	Lyophilized
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	10 mM Ammonium Bicarbonate (NH ₄ HCO ₃), pH 7.4.
Preservative	None
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Apolipoprotein A4;Apolipoprotein antigen;Apolipoprotein protein