



Recombinant human apolipoprotein E3 (>90%) (DAGA-961)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Apolipoprotein A1 is a protein that in humans is encoded by the APOA1 gene. It has a specific role in lipid metabolism. The text in a recent report suggested that APOA7 mRNA is regulated by endogenously expressed antisense RNA.
Purity	>90% by SDS-PAGE and HPLC
Conjugate	TBD
Applications	Specific methodologies have not been tested using this product
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	20 mM Sodium Phosphate, pH 7.8 containing 0.5 mM DTT
Preservative	None
Storage	Store lyophilized product at -20°C. After reconstitution, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Apolipoprotein A7;Apolipoprotein antigen;Apolipoprotein protein