



Human fibrinogen (DAGA-965)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Fibrinogen (factor I) is a glycoprotein in vertebrates that helps in the formation of blood clots. It consists of a linear array of three nodules held together by a very thin thread which is estimated to have a diameter between 8 and 15 Angstrom (Å).
Conjugate	Unconjugated
Applications	Specific methodologies have not been tested using this product
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Buffer	10 mM Sodium Citrate, 10 mM Sodium Phosphate, pH 7.3.
Preservative	None

BACKGROUND

Keywords	Fibrinogen;factor I
-----------------	---------------------