



B. anthracis Edema Factor (DAGB168)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	B. anthracis
Purity	>95%
Conjugate	Unconjugated
Applications	ELISA
Molecular Weight	89kDa
Reconstitution	When reconstituted with 0.1 ml of water, the concentration of buffer is 5mM HEPES, 50mM NaCl, pH7.5.
Endotoxin	The endotoxin content, determined using a kinetic chromagenic LAL assay is 0.1 EU/mg.
Format	Lyophilized power, sealed under vacuum
Size	0.1 mg
Preservative	None
Storage	Store at 2-8°C prior to reconstitution.