



Eastern Equine Encephalitis Virus E3E2 Protein [His] (DAGB175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Eastern Equine Encephalitis Virus
Purity	>90% determined by SDS-PAGE
Conjugate	His
Applications	WB
Molecular Weight	55 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg, 1 mg
Buffer	20mM Tris-HCl buffer, pH 8.0, 300 mM NaCl, 7M Urea
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.
Ship	Dry ice