



# Western Equine Encephalitis Virus E2 GP [His] (DAGB176)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Species</b>	Western Equine Encephalitis Virus
<b>Conjugate</b>	His
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Concentration</b>	0.615 mg/mL
<b>Size</b>	100 µg
<b>Buffer</b>	PBS (supplemented with glycerol, arginine and glutamic acid)
<b>Preservative</b>	None
<b>Storage</b>	2-3 weeks at -20°C, -80°C long term