



Bundibugyo Glycoprotein (minus the Transmembrane Region) [His] (DAGB178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bundibugyo Ebolavirus
Conjugate	His
Applications	WB
Molecular Weight	~73 kDa including the His-tag, without glycosylation
Format	Liquid
Concentration	1.59 mg/mL
Size	100 µg
Buffer	PBS (supplemented with 10% glycerol, arginine and glutamic acid)
Preservative	None
Storage	2-3 weeks at -20°C, -80°C long term