



Angola marburgvirus Glycoprotein (minus the Transmembrane Region) [His] (DAGB185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Angola marburgvirus
Conjugate	His
Applications	ELISA, WB
Molecular Weight	~60 kDa including the His-tag, without glycosylation
Format	Liquid
Concentration	0.55 mg/mL
Size	100 µg
Buffer	PBS (supplemented with glycerol, arginine and glutamic acid)
Preservative	None
Storage	2-3 weeks at -20°C, -80°C long term