



Recombinant IAV H10N8 (A/Jiangxi-Donghu/346/2013) Hemagglutinin (DAGB189)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	IAV H10N8
Conjugate	Unconjugated
Applications	WB
Molecular Weight	~62 kDa, without glycosylation
Format	Liquid
Concentration	0.418 mg/mL
Size	50 µg
Buffer	20 mM Tris pH 7.4, 150 mM NaCl and 0.01% non-ionic surfactant
Preservative	None
Storage	2-3 weeks at -20°C, -80°C long term

BACKGROUND

Keywords Influenzavirus A; Influenza A virus; Influenza A virus H10N8 HA; H10N8; HA; IAV H10N8 HA