



Salinomycin [HRP] (DAGB312)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---|
| Specificity | Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development. |
| Species | N/A |
| Conjugate | HRP |
| Applications | IA |
| Format | The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution |
| Size | 0.5 ml |
| Preservative | None |
| Storage | 2 - 8°C for up to 3 months / -20°C for long term storage |

BACKGROUND

| | |
|---------------------|--|
| Introduction | Salinomycin is an antibacterial and coccidiostat ionophore therapeutic drug. |
|---------------------|--|