



Salinomycin [OVA] (DAGB312B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	0.01M pH7.4 PBS
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Salinomycin and its derivatives exhibit high antimicrobial activity against Gram-positive bacteria, including the most problematic bacteria strains such as methicillin-resistant <i>Staphylococcus aureus</i> and methicillin-resistant <i>Staphylococcus epidermidis</i> , and <i>Mycobacterium tuberculosis</i> . Salinomycin is inactive against fungi such as <i>Candida</i> and Gram-negative bacteria.
Keywords	Salinomycin; Coccidiostat