



Maduramicin [HRP] (DAGB313)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Size	0.5 ml
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage

BACKGROUND

Introduction	Maduramicin is a monoglycoside polyether derived from the fungus <i>Actinomadura rubra</i> . It is classified as an ionophor and predominantly used as an anti-coccidial agent, which is mixed with poultry feed in the ratio of 6ppm.
---------------------	--