



Desmethylnaproxen [HRP] (DAGB329)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxide and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage

BACKGROUND

Introduction	Desmethylnaproxen is an anti-inflammatory agent with analgesic and antipyretic properties. Both the acid and its sodium salt are used in the treatment of rheumatoid arthritis and other rheumatic or musculoskeletal disorders, dysmenorrhea, and acute gout.
---------------------	--