



Dicloxacillin [HRP] (DAGB335)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage

BACKGROUND

Introduction	Dicloxacillin is a narrow-spectrum beta-lactam antibiotic of the penicillin class. It is used to treat infections caused by susceptible Gram-positive bacteria. It is active against beta-lactamase-producing organisms such as Staphylococcus aureus, which would otherwise be resistant to most penicillins. It is very similar to flucloxacillin and these two agents are considered interchangeable.
---------------------	--