



Norfentanyl [HRP] (DAGB340)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Size	0.5 ml
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage

BACKGROUND

Introduction	Norfentanyl is a metabolite of Fentanyl. It is a potent narcotic analgesic, abuse of which leads to habituation or addiction. It is primarily a mu-opioid agonist. Norfentanyl is also used as adjunct to general anaesthetics, and as an anaesthetic for induction and maintenance.
---------------------	--