



Tetrahydrocorticosterone [HRP] (DAGB352)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxidase and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage

BACKGROUND

Introduction	3 α ,5 α -Tetrahydrocorticosterone (3 α ,5 α -THB), or simply tetrahydrocorticosterone (THB or THCC), is an endogenous glucocorticoid hormone.
---------------------	--