



3-Acetyl-DON [HRP] (DAGB356)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Each conjugate comprises antigen covalently bound to horseradish peroxide and is suitable as a tracer in immunoassay development.
Species	N/A
Conjugate	HRP
Applications	IA
Format	The conjugate is supplied as a concentrate. Dilute as required and use working strength conjugate immediately after dilution
Preservative	None
Storage	2 - 8°C for up to 3 months / -20°C for long term storage