



## NIP [BSA] (DAGB485)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                                                                                       |
|--------------|---------------------------------------------------------------------------------------|
| Species      | N/A                                                                                   |
| Conjugate    | BSA                                                                                   |
| Applications | ELISA                                                                                 |
| Format       | Liquid                                                                                |
| Size         | 10 mg, 100 mg                                                                         |
| Buffer       | PBS buffer                                                                            |
| Preservative | None                                                                                  |
| Storage      | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| Ship         | Ambient                                                                               |

### BACKGROUND

|              |                                                                                                         |
|--------------|---------------------------------------------------------------------------------------------------------|
| Introduction | 4-Hydroxy-3-iodo-5-nitrophenylacetyl hapten is conjugated to BSA protein through lysine by amide bonds. |
|--------------|---------------------------------------------------------------------------------------------------------|