



## NIP [BSA-Biotin] (DAGB487)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	N/A
Conjugate	BSA-Biotin
Applications	ELISA
Format	Liquid
Size	5 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

### BACKGROUND

Introduction	Biotin O-Succinimide is conjugated to 4-Hydroxy-3-iodo-5-nitrophenylacetyl-Bovine Serum Albumin.
--------------	--