



NP [BSA], Ratio 1-9 (DAGB493)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA
Applications	ELISA
Reconstitution	Reconstitute with deionized sterile water
Format	Lyophilized
Size	10 mg, 100 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	4-Hydroxy-3-nitrophenylacetyl hapten is conjugated to BSA (Bovine Serum Albumin) protein through lysine by amide bonds.
--------------	---