



NP [BSA-Fluorescein] (DAGB494)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA-Fluorescein
Applications	ELISA
Format	Liquid
Size	5 mg, 10 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction	Fluorescein is conjugated to 4-Hydroxy-3-nitrophenylacetyl-BSA through lysine by amide bonds.
--------------	---