



NP [BSA-Biotin] (DAGB495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA-Biotin
Applications	ELISA
Format	Liquid
Size	5 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction	Biotin is conjugated to NP-BSA.
--------------	---------------------------------