



## NP [CGG], Ratio 10-19 (DAGB497)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	N/A
Conjugate	CGG
Applications	ELISA
Reconstitution	Reconstitute with deionized sterile water
Format	Lyophilized
Size	1 mg, 5 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

### BACKGROUND

Introduction	4-Hydroxy-3-nitrophenylacetyl hapten is conjugated to CGG lysine through amide bonds.
--------------	---