



NP [CGG], Ratio 30-39 (DAGB500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	CGG
Applications	ELISA
Format	Liquid
Size	1 mg, 5 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction	4-Hydroxy-3-nitrophenylacetyl hapten is conjugated to CGG lysine through amide bonds.
--------------	---