



## NP [CGG-Fluorescein] (DAGB501)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	N/A
Conjugate	CGG-Fluorescein
Applications	ELISA
Format	Liquid
Size	1 mg, 5 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

### BACKGROUND

Introduction	Fluorescein is conjugated to NP-CGG lysine through amide bonds.
--------------	---