



NP [Dextran] (DAGB502)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	Dextran
Applications	ELISA
Format	Liquid
Size	10 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction	NP (4-Hydroxy-3-nitrophenylacetic) hapten conjugated to Amino DEXTRAN. Conjugation is given as NPs per molecule of DEXTRAN and is determined by NP absorbance at 430 nm, assuming DEXTRAN MW of 500,000.
--------------	--