



NP [Hexyl-Fluorescein] (DAGB505)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	Hexyl-Fluorescein
Applications	ELISA
Molecular Weight	653.64
Format	yellow powder
Size	1 mg, 5 mg
Buffer	Dissolve in DMF or DMSO. Store solution in aliquots at -20°C. Protect from light.
Preservative	None
Storage	2 to 8 °C
Ship	Ambient

BACKGROUND

Introduction 4-Hydroxy-3-nitrophenylacetyl diaminohexane is conjugated to fluorescein via amide linkage.