



NP [LPS] (DAGB509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	LPS
Applications	ELISA
Format	Liquid
Size	1 mg, 5 mg
Buffer	PBS buffer
Preservative	None
Storage	2 to 8 °C
Ship	Ambient

BACKGROUND

Introduction	4-Hydroxy-3-nitrophenylacetyl hapten is conjugated to LPS (Lipopolysaccharide) through lysine by amide bonds.
--------------	---