



NP [PE] (DAGB511)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	PE
Applications	ELISA
Format	Liquid
Size	1 mg, 5 mg
Buffer	PBS, pH7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	4-Hydroxy-3-nitrophenylacetyl hapten is conjugated to PE (Phycoerythrin) protein through lysine by amide bonds.
Keywords	4-Hydroxy-3-nitrophenyl; NP; PE; Phycoerythrin