



TNP [BSA] (DAGB523)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	BSA
Applications	ELISA
Format	Liquid
Size	10 mg, 100 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction	2,4,6-Trinitrophenyl hapten is conjugated to Bovine Serum Albumin protein through lysine by amide bonds.
--------------	--