



TNP [KLH] (DAGB526)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	KLH
Applications	ELISA
Format	Liquid
Size	5 mg, 25 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction 2,4,6, Trinitrophenyl hapten is conjugated to KLH protein through lysine by amide bonds.