



TNP [KLH-Biotin] (DAGB527)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	N/A
Conjugate	KLH-Biotin
Applications	ELISA
Format	Liquid
Size	5 mg
Buffer	PBS buffer
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Ship	Ambient

BACKGROUND

Introduction	Biotin OSu is conjugated to TNP-KLH. This product can be used for ligand-receptor interactions in cell cytometry.
--------------	---