



Recombinant Hepatitis B virus Surface Antigen S-protein (adr) (DAGC049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------|---|
| Purity | > 95% as determined by SDS-PAGE |
| Conjugate | Unconjugated |
| Applications | ELISA, WB |
| Format | Liquid |
| Size | 100 µg |
| Buffer | 10mM sodium phosphate, 140mM NaCl, pH 6.4 |
| Preservative | 0.1% sodium azide |
| Storage | Store at 4°C |

BACKGROUND

| | |
|-----------------|---|
| Keywords | HBsAg S-protein; HBsAg; Small HBV surface protein |
|-----------------|---|