



Recombinant Human Proinsulin Protein (a.a. 25-110) [His] (DAGC108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	≥95% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, FCM
Molecular Weight	10.5 kDa
Reconstitution	Reconstitute at 100 µg/mL in sterile PBS
Format	Lyophilized
Size	50 µg
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Proinsulin; Insulin; C-peptide
----------	--------------------------------